


LevyLab Instrument Framework


Patrick Irvin

June 29, 2020

Experiment Overview



Front Panel (FP) Valve Controller




Dilution
Refrigerator
(T=20mK)



Temperature Controller (TC)



LabVIEW Instrument (FP)




LabVIEW Instrument
(PS)



Superconducting
Magnet Power Supply (PS)

Superconducting
Magnet



sample




LabVIEW Instrument (TC)




Data Acquisition (DAQ)

Software Frameworks and Design Patterns


LabVIEW



LabVIEW...




LabVIEW...



Software Frameworks Can Help Us!


- **JKI State Machine:** A powerful state machine design template
- **Object Oriented Design:** An approach to designing a system of interacting objects
- **State Machine Objects:** Turn JKI State Machines into asynchronous objects!
- **LevyLab Instrument.lvclass:** A framework using JKI State Machines, SMO, that solves specific needs of LevyLab Research




Software Frameworks and Design Patterns

JKI State Machine

JKI State Machine




JKI State Machine




JKI State Machine

Event example:




JKI State Machine

State example:



JKI State Machine Explorer



JKI State Machine Best Practices

JKI State Machine Best Practices

<http://blog.jki.net/products/state-machine/jki-state-machine-best-practices>


1. [Don't hide your state strings in subVIs](#)
2. [Don't add code and logic inside the Event Structure](#)
3. [Keep the Original Size \(i.e. don't grow the structures\)](#)
4. [Use macros instead of “chaining” together sequential states](#)
5. [Left-justify State Strings instead of Right-justify](#)

Software Frameworks and Design Patterns

Object Oriented Programming with LabVIEW

Object Oriented Programming LabVIEW


- [NI OO Course](#)
- <https://www.bloomy.com/support/blog/object-oriented-labview-inheritance-part-1-3-part-series>
- <https://www.bloomy.com/support/blog/object-oriented-labview-encapsulation-part-2-3-part-series>
- <https://www.bloomy.com/support/blog/object-oriented-labview-polymorphism-part-3-3-part-series>




Software Frameworks and Design Patterns

State Machine Objects


State Machine Objects




State Machine Objects



State Machine Objects: Dynamic Events




Awesome! How do I make an SMO?



1. Open the SMO Editor





2. Click "New SMO"




3. Choose SMO Type

SMO Intro 1




SMO Intro 2a




Startup:

Test Launcher




Process.vi




During Execution:


Calling PublicMethod from TestLauncher triggers event and sends data to process.vi




Calling PublicEvent from Process.vi triggers event and sends data to TestLauncher



SMO Intro 2b




Use SMO Editor to create Public Events and Public Methods




SMO Intro 3

Change your icons to make them useful!



Use glyphs provided by LabVIEW to make your life easy!



SMO Resources

- [JKI State Machine Objects](#) on GitHub
 - [Getting Started with JKI's State Machine Objects Framework](#)
- [Patrick's Short Introduction to JKI SMO](#)
- LevyLab [SMO Tutorials](#)


Software Frameworks and Design Patterns

Instrument.lvclass


Instrument Framework

- <https://github.com/levylabpitt/Instrument-Framework>
- An Instrument developed using this framework will have access to the following capabilities:
 - Methods for reading and writing configuration settings.
 - A hardware abstraction layer (HAL) for interfacing with hardware.
 - Data logging: poll and log instrument configuration and settings to DSC database
 - In development: logging to Amazon AWS
 - Defines an API to allow external programs to control compiled instances of the application and remotely across a network
 - Currently implemented using [National Instruments' Simple Messaging Library \(STM\)](#), however a cross-platform protocol is being developed using [OMQ](#).
 - User interface framework.
 - In development: an embedded subpanel UI framework ([MAUI](#))
- The LevyLab Instrument Framework makes extensive use of [JKI State Machines](#) and [JKI State Machine Objects](#).


Instrument Framework Diagram



Instrument.lvclass



Instrument.lvclass




Instrument Process (SMO): The main process that runs behind the scenes

- Starts OTHER SMOs (Loggers (DSC/Error), Remote Control Server)
- Periodically instrument state to DSC
- Handles messages from Remote Server
- Provides configuration framework


Instrument.lvclass

Handle Command.vi

- Overridable method that knows what to do with an API call




Instrument.lvclass




Remote Control (SMO)

- Knows how to communicate over network connection (TCP, 0MQ, etc)
- Manages remote connections
- Brokers commands between instrument and client


Instrument.lvclass



Instrument.lvclass: UML Diagram



Instrument.lvclass: UML Diagram (Responsibilities)



Create your own Instrument

Instrument.lvclass Prerequisites

- Packages:
 - JKI State Machine Objects (SMO)
 - JKI SMO Editor
 - LevyLab Instrument Template (Base with UI)



Create Instrument. *NewInstrument.lvclass*

1. Setup repository in Bitbucket (\levylab_inst_NewInstrument)
2. Setup codebase using the SMO Editor
3. Modify overridable methods
4. Inside AppLauncher.vi:
 1. Your Instrument must have an SMO Name (HWName_S/N)
 2. ON HOUSESOFTHEHOLY: Create a Process in the DSC that matches the SMO Name above (example: SR7270_09337296)
 3. Change class object to your NewInstrument
 4. ~~#TODO: "CreateSMOName.vi" should be part of instrument class. child must override~~
 5. ~~#TODO: ensure children must override CERTAIN protected methods~~
5. getAll
 1. getAll data must correspond to variables inside the DSC
 2. ON HOUSESOFTHEHOLY: create variables with logical names that match the getAll datatype cluster
6. Make a UI

1. Setup repository in Bitbucket (`\levylab_inst_NewInstrument`)

1. Option 1: If this will be a new instrument, create a new repository on Bitbucket and clone to a working directory on your computer (ie, `c:\user\levylab_inst_NewInstrument`)
2. Option 2: If you are refactoring an existing instrument, you can make a new branch and work from there.

2. Setup codebase using the SMO Editor



3. Modify overridable methods

1. Modify these VIs:

1. Create Instrument SMO Name.vi
2. Open Instrument.vi
3. Close Instrument.vi
4. Handle Command.vi
5. getAll.vi

2. Modify these Typedefs:

1. *NewInstrument.Commands--Enum.ctl*
2. *NewInstrument.Configuration--Cluster.ctl*
3. *NewInstrument.getAll--Cluster.ctl*


3. Optional

1. Configure Instrument.vi (and subVIs)

4. Other (do not need to modify)

1. Get Public API


1.1 Create Instrument SMO Name.vi



Your SMO name needs to match an existing process in the DSC


1.2 Open Instrument.vi

- Called by Instrument.lvclass:Process.vi's state Remote: Init (called by Macro: Initialize)
- Add code to open connection to hardware
- Example:



1.4 Handle Command.vi

- Modify code to respond to command strings
- Could talk to hardware directly or send Public Method to another SMO
- Example:



Typedefs

1. `NewInstrument.Commands--Enum.ctl`


an enum containing a list of supported API commands



2. `NewInstrument.Configuration--Cluster.ctl`


any special configuration settings your instrument needs.

Hardware address?



3. `NewInstrument.getAll--Cluster.ctl`

this is a cluster that will be written to the DSC. The name of each element needs to match an existing trace in the DSC




PPMS2


- Instrument.PPMS2.lvclass already exists... What is the plan?
 - Copied Instrument.Example and *.Example.UI merge code....
- Make existing *.PPMS2.lvclass inherit from Instrument.lvclass >v1.1.0.9 **check**
- Make inherit from Instrument.Cryostat.lvclass **check**
- Make new overrides (set Field, etc) **check**
 - Get Angle, Get Field, Get Helium Level, Get Temperature, Set Angle, Set Field, Set Temperature
- Override “~~Initialize~~ Open Instrument”
 - Copy code from “PPMS2: Initialize” state (?) **check**
 - Add Instrument Ref to Instrument.PPMS2 class private data **check**
- Override “Get Public API” **check**
- Override “Handle Command” **check**
 - getAll: code from “PPMS2: Get Status” state **check**
 - (make override method “getAll”) **check**
- Update APIs to use remote client **check**
- Override “Create Instrument SMO Name.vi”
 - Return name “PPMS 2” **check**

Overrides


Open Instrument.vi



getAll.vi




Get Public API.vi




Create Instrument SMO Name.vi

API & Override “Handle Command.vi” (Examples)

Set Magnet Field.vi (API)




Handle Command.vi: Set Magnet



Non-Standard Stuff...



- Instrument.PPMS2: PPMS Monitor and Control.vi
 - Keep this one, but make code match Instrument.Example:AppLauncher.vi
- Instrument.PPMS2: Process.vi
 - Ultimately I'm not sure if this one should be doing anything. Move UI elements to →
- Instrument.PPMS2 UI.lvclass
 - Create this class. This will contain the UI
 - This class already exists as Instrument.Example UI.lvclass. Just rename that stuff

AnyInstrument.lvclass ↔ Instrument.lvclass




AnyInstrument.lvclass:process.vi

- Override “Handle Remote Command” and “Get Public API” to list of supported commands
- Commands are processed in “Handle Remote Command”
- If Handle Remote Command needs to communicate with AnyInstrument’s process it uses methods provided by AnyInstrument



1. Instrument process receives event from **RemoteControl** and calls **Handle Remote Command**
2. **Handle Remote Command** decides what to do with command. In this example it called **AnyInstrument:setMagnet**
3. **AnyInstrument:setMagnet** fires event
4. **setMagnet** event is received by **AnyInstrument**.
5. **AnyInstrument** does its work and *responds using notifier** provided by **setMagnet***
6. **setMagnet** receives notifier from **AnyInstrument**
7. **Instrument** calls **RemoteControl:SendMessageToProcess** to send processed data



1. **Handle Remote Command*** is called by **Instrument’s process**.

Is an overridable method by child class.

7. Any responses are forwarded to **RemoteControl** by calling **RemoteControl.CommandResponse**


*Handle Remote Command is protected.

Can also have “Script Remote Command” that is a public method

**This programming structure provides an easy way for instrument to respond to the API when it is finished




Instrument.lvclass: EnumerateStaticDependencies.vi

EnumerateStaticDependencies







Make `RemoteControl.lvclass` a dependency of `Instrument`
Call `write` commands and port here to initialize `RemoteControl`


`EnumerateStaticDependencies` called here



Instrument.lvclass: Handle Remote Command



Instrument.lvclass ↔ RemoteControl.lvclass




Instrument.lvclass:process.vi

- Uses **RemoteControl** as dependency (declared by `enumerateStaticDependencies.vi` (overridden from SMO))
- Processes commands received from **RemoteControl.SendMessageFromProcess**
- List of supported commands maintained by overridable method "Get Public API"
 - children of `Instrument` override with their own list
- Commands are processed by overridable method "Handle Remote Command"
- Responds by calling **RemoteControl.SendMessageToProcess**

RemoteControl.SendMessageFromProcess

RemoteControl.SendMessageToProcess



RemoteControl.lvclass:process.vi

- Acts as command server between external processes (clients) and owner process
- Should be protocol agnostic (should work for raw TCP, STM, OMQ by choice of child class)


Methods are owned by **RemoteControl**.

RemoteControl.lvclass \longleftrightarrow RemoteControl (.STM) or (.0MQ)

Instrument

RemoteControl.SendMessageFromProcess

RemoteControl.SendMessageToProcess



Read/Write

Write/Read


STM Client
---or---
0MQ Client

RemoteControl.SendMessageFromProcess

RemoteControl.SendMessageToProcess


RemoteClient

RemoteControl Diagram (11/29/2018)





RemoteControl Description (11/29/2018)


1. “Instrument” starts a RemoteControl (RC) Connection Monitor SMO
2. Connection Monitor listens for clients to connect
3. “Client” starts a RC Client SMO
4. Connection Monitors receives the connection from Client and starts a RC Server SMO
5. Both RC Server and Client poll Read Message according to the timeout
6. Client sends a message by generating a private user event. Client process catches event and calls Client’s Send Message to send message over TCP to Server’s Read Message (5).
7. Server generates public event that is caught by Instrument Process.
8. Instrument processes message and generates private user event to talk back to Client through same call chain described in 6&7.
20. (*how does Instrument know about Server or vice versa?*)




RemoteControl Diagram (6/26/2019)



RemoteControl Tests




1. Start Connection Monitor
2. Client send 100 messages of 100 B size and reports average number of messages/second, time/message, and throughput (MB/s)
3. Increase message size and repeat up to ~10MB



Discussion with Javier 6/27/2019


- RemoteControl Issue:
 - Server process would only stop half of the time
 - If Client closes the connection, Server receives error 66, Server process calls Macro: RC Exit, Macro: RC Exit calls SelfTerminate.
 - SelfTerminate was not being called every other time, because Macro: RC Exit was being called before a selector controlled by “Core.Close Panel on Exit?”




Discussion with Javier 6/27/2019

- Build issue
 - \Templates\ExampleInstrument\ and \Templates\ExampleInstrumentUI\ have some cross-linking. They need to be build in one package
 - ?? – can ExampleInstrumentUI.lvclass be named ExampleInstrument.UI.lvclass? Maybe this would fix the issue with needed to change the class wires after the class is created from the template.

RemoteControl (RC) Diagram




UI



Instrument.Cryostat.lvclass: API Description

- Magnet
 - Set Magnet Field[Target Field, Target Rate, Mode, Axis]
 - Target Field: dbl (Teslas)
 - Target Rate: dbl (Teslas/minute)
 - Mode: enum: 0: Persistent, 1: Driven
 - Axis: enum: 0: Z (*Default*), 1: X, 2: Y
 - Get Magnet Field[Axis] → **Field, Status**
 - Axis: enum: 0: Z (*Default*), 1: X, 2: Y
 - **Field: dbl (Teslas)**
 - **Status: variant**
 - Rotator
 - Set Angle[Angle, Axis]
 - Angle: dbl (°)
 - Axis: enum: 0, 1
 - Get Angle[Axis] → **Angle**
 - Axis: enum: 0, 1
 - **Angle: dbl (°)**
 - Temperature
 - Set Temperature[Temperature, Rate, Channel]
 - Temperature: dbl (K)
 - Rate: dbl (K/min)
 - Channel: int
 - Get Temperature[Channel] → **Temperature, Status**
 - Channel: int
 - **Temperature: dbl (K)**
 - **Status: variant**
 - Level
 - Get Helium Level → **Level**
 - **Level: dbl (%)**
 - Get Nitrogen Level → **Level**
 - **Level: dbl (%)**
- Command[argument, [optional argument]] → **Response**
 - An Instrument does not need to override all methods
 - Example 1: PPMS/PPMS2 Model6000 handles all of these commands, so override them all in `Instrument.Cryostat.PPMS2`
 - Example 2: MNK system uses separate instruments for each:
 - `Instrument.Cryostat.Attocube` overrides:
 - Rotator
 - `Instrument.Cryostat.OxfordITC` overrides:
 - Temperature
 - `Instrument.Cryostat.OxfordIPS` overrides:
 - Level, Temperature, Magnet
 - `Instrument.Cryostat.LeidenTC` overrides:
 - Temperature
 - `SweepControl.vi` allows different instruments for each type (Magnet, Temperature, etc)



- There was a problem building PPMS2 Instrument.
- Error Details:
- Code: 1357
- Source: Error 1357 occurred at Invoke Node in Save VIs to Destination __ogb.vi->Build Application __ogb.vi->Build Application from Build File __ogb.vi->Build Application from Build File __ogb.vi.ProxyCaller
- Possible reason(s):
 - LabVIEW: A LabVIEW file from that path already exists in memory, or exists within a project library already in memory.
- Error Conditions:
- Private Save
- VI Name: "7F5BF585CD3CEE1343C47873D09E4DD8Logger.DSC.lvclass:PrivateEvents--Cluster.ctl"
- VI Path: "C:\Users\Patrick\AppData\Local\Temp\vpbttmp\src\LabVIEW\user.lib\LevyLab_PPMS2Instrument_internal_deps\7F5BF585CD3CEE1343C47873D09E4DD8PrivateEvents--Cluster.ctl"
- Method Name: Save:Instrument (Private)

VI PPMS2 Instrument - [PPMS2 Instrument.vipb] - VI Package Builder

File Edit Package Tools Window Help

Basic

- Build Information
- Display Information
- Palettes
- Destinations
- Source File Settings

Advanced

- Package Dependencies
- Licensing & Activation
- Install Requirements
- Incompatible Packages
- Custom Actions
- Pre/Post Build Actions
- Package Filename

Package Dependencies

Package Name	Version	Description
JKI State Machine Objects (SMO)	>=1.3.0.56	Source VIs call VIs in this package
OpenG Application Control Library	>=4.1.0.7	Source VIs call VIs in this package
OpenG Time Library	>=4.0.1.3	Source VIs call VIs in this package
Instruments	>=1.2.4.37	User Added

VI 2016

+ - X

File Edit View Package Tools Window Help


VI NEW VI FOLDER VI DOCUMENT VI CHECK VI DOWNLOAD

2016

Name / Version Repository Company

JKI State Machine	2018.0.7.45	JKI Package Network	JKI
JKI State Machine Objects (SMO)	1.3.0.56	JKI Package Network	JKI
jkI_rsc_toolkits_palette	1.1-1	JKI Package Network	JKI Software
OpenG Application Control Library	4.1.0.7	JKI Package Network	OpenG.org
OpenG Error Library	4.2.0.23	JKI Package Network	OpenG.org
OpenG LabVIEW Data Library	4.2.0.21	JKI Package Network	LAVA
OpenG Time Library	4.0.1.3	JKI Package Network	National Instruments
Instruments			


Ready ...



Instrument.VNA.lvclass: API Description

- The screenshot displays two VI front panels in LabVIEW:

 - Agilent P9374A VI:** This VI is organized into several sub-VIs and a main configuration area. Sub-VIs include "Mag", "S11", "S21", "S12", "S22", "Phase", and "Data". The main configuration area contains controls for "power", "start/center", "stop/span", "IF bandwidth", "average", "average factor", "saver's name", and "graphical display". It also includes a "VISA resource name" input, an "error in (no error)" output, and a "VISA resource name" output.
 - Agilent HP8753D VI:** This VI is presented as a grid of function icons under the heading "VIs for HP8753D (Dutt's)". The grid is organized into columns corresponding to the LabVIEW menu bar: Initialize, Configuration, Action/Status, Data, Utility, and Close. Each column contains multiple sub-functions, such as "Initialize", "Configuration", "Action/Status", "Data", "Utility", and "Close", each with its own specific set of icons.



Instrument.VNA.lvclass: API Description

Step 1: Set VNA parameters
Directly give commands to VNA

- Initialize VNA
- Set Power
- Set measurement
 - S11, S21, S12, S22 (choose)
- Set format
 - Log Mag, Phase (choose)
- Set sweep
 - Start/stop or center/span (T/F)
 - Start/center frequency
 - Stop/span frequency
 - Number of points
 - Sweep time
- Set average
 - Average or not(T/F)
 - Average factor
 - IF bandwidth
- Collect data setting
 - Save to file (T/F)
 - Data format
 - Folder, file, title, note...

Step 2: Sweep control and take VNA data every sweep

sweep bias
or sweep Vbg
or sweep B field
or sweep temperature
or ...
by Sweep control


Collect a sweep data from VNA
(smallest loop)

- Frequency
- S parameters
- Log mag, phase


Instrument.VoltageUpdate.lvclass

- VoltageUpdate private data:
 - ~~+ gain in “configuration”~~
 - Add a configuration window
 - ~~Voltage → DAQ Voltage & Amplified Voltage~~
- Constants
 - SMO Name-Constant.vi (overridable)
 - ~~SupportedCommands Constant.vi (overridable)~~
- ~~Configure API Methods~~
- ~~VoltageUpdate UI~~
- DSC Process
 - Option 1: PPMS 2 VoltageUpdate
 - Option 2: VoltageUpdate_ComputerName (e.g. VoltageUpdate_BLACKPEARL)
 - Option 3: VoltageUpdate_S/N
- Instrument.lvclass:
 - Add accessor DSC process
 - Joe Albro: VI tags for scripting Instrument/Instrument UI.vi
 - Configuration stuff is a mess


Experiment 1: Hardware Overview



Front Panel (FP) Valve Controller



LabVIEW Instrument (FP)



LabVIEW Instrument (PS)



Superconducting Magnet Power Supply (PS)

Dilution Refrigerator (T=20mK)




Temperature Controller (TC)



LabVIEW Instrument (TC)




sample




Data Acquisition (DAQ)

Experiment 2: LabVIEW Instrument

- A LabVIEW Instrument has the following responsibilities:
 - Know how to communicate with a piece of hardware (drivers)
 - Poll the instrument for its settings and log to a database
 - Open an API for external programs (e.g through cross-platform protocol such as 0MQ)
 - Provide a UI (optional but probably desirable)



Experiment 3: An “Experiment”





1. Configure wiring and voltage and current amplifiers

2. Define global settings such as sample name and description.

3. Define experiment

Initialize some stuff
For $B = 0$ to B_{\max} :
 Set magnet = B
 For $V_{\text{drive}} = 0$ to $V_{\text{drive},\max}$:
 For $i = 0$ to N
 Get current_i
 Get voltage_i
 Write current and voltage to File_B.itzx
 Exit



Experiment 4: ITX (Igor Text File)* Description 1

```
1 IGOR
2 X// Date: 1/12/2017 8:40 PM
3 X// sweep B: 3.5T to 9.5T, 0.06 T/min
4 X// large axis: 147, small axis: 257 (20 deg)
5 X// sweep Vsg: 0mV to 120mV 250uV step
6 X// Vbg=0.3V
7 X//
8 X// tc=0.3s
9 X// order=4
10 X// 60 Hz filter OFF
11 X//
12 X// lockin:
13 X// 0: 1, voltage
14 X// 1: 2, current, 100uV AC 240uV DC, 13 Hz
15 X// 2: 3, voltage
16 X// 3: 4, voltage
17 X// 4: 5, current
18 X// 5: 7, voltage
19 X//
20 X// KH gain=100x, shunt resistor=50kOhm
21
```

- First line is always “IGOR”

- Metadata begins with X//
- This section should describe the experiment that was done
- This one describes that 1. B will be swept from 3.5 Tesla to 9.5 Tesla at a rate of 0.06 Tesla/minute. 2. Vsg will be swept from 0 mV to 120 mV in 250 uV steps.
- Each file will be one Vsg sweep. The folder will contain data between B = 3.5 T and 9.5 T. And as many Vsg sweeps that could be taken while B was changing.

- This section describes how the lockin was configured.
- The lockin measures current or voltage on the sample in response to some source voltage stimulus.
- There are typically 8 lockin channels (not all need to be used) and there are X and Y for each, which are in-phase and out-of-phase, respectively, with the source voltage.
- The information here means:
 - Lockin0 (connected to electrode 1) measures a voltage
 - Lockin1 (connected to electrode 2) sources a current using a 100 uV source voltage at 13 Hz
 - ...
 - Lockin5 (connected to electrode 7) measures a voltage

Experiment 4: ITX Description 2

```
21  
22 WAVES/D/N=(481) 'time,SA02836H.20170110.000002'  
23 BEGIN  
24 3.567116343758267E+9  
25 3.567116353762839E+9  
26 3.567116353912848E+9  
27 3.567116354064857E+9  
28 3.567116354215865E+9  
29 3.567116354440878E+9  
30 3.567116354590887E+9  
31 3.567116354818900E+9
```

- Here starts a wave of size 481. This one is the time
- (seconds since epoch time (1904 I think))

```
497 3.567116453997572E+9  
498 3.567116454222585E+9  
499 3.567116454298590E+9  
500 3.567116454448598E+9  
501 3.567116454599607E+9  
502 3.567116454824620E+9  
503 3.567116454986629E+9  
504 3.567116455136638E+9  
505 END
```

- Sample code: SANNNNNL.YYYYMMDD
 - SA: Sample
 - NNNN: 5 digit sample number
 - L: one of 16 “canvases” on a particular sample
 - YYYY: Year
 - MM: Month
 - DD: Day

- End of the wave

```
1477 WAVES/D/N=(1) 'B.SA02836H.20170110.000002'  
1478 BEGIN  
1479 3.801700000000000E+0  
1480 END
```


- Wave of size 1 showing the B field for THIS file

```
1481  
1482 WAVES/D/N=(481) 'X0.SA02836H.20170110.000002'  
1483 BEGIN  
1484 100.586757640980660E-6  
1485 100.562881961642580E-6  
1486 100.554882945805970E-6  
1487 100.548795103440390E-6  
1488 100.541729931379820E-6  
1489 100.540689133795410E-6  
1490 100.545901069673990E-6  
1491 100.563244064805010E-6  
1492 100.589700051043400E-6
```

- Wave for In-phase (X) Lockin0

Experiment 5: Data Processing 1

1. Look at experiment notebook for sample SA02836H.20170110. There should be a description & diagram of the device:



2. We want to plot the Four terminal resistance of the exciting section. So we need a current and voltage on each side of the device.

2a. The metadata tells us electrode 2 and 5 were current source and drain. Then pick a voltage on each side, e.g. 3 and 4, 1 and 7, 3 and 7, or 1 and 4.

2b. Now go back to the metadata again and decipher which lockins those correspond to*. Let's choose electrodes 1(V), 2(I), 4(V), 5(I), which are lockins 0(V), 1(I), 3(V), 4(I).


3. The 1st dependent variable is V_{sg} . The 2nd is B

4. Read the waves X0.* , X1.* , X3.* , X4.* , Vsg.* , B.* from each file. Calculate $R_{4T} = \frac{V_1 - V_2}{I} = \frac{X3 - X0}{X4}$. You now have R_{4T} vs V_{sg} for each file, each of which is one value of B .


*I will not be offended if you roll your eyes at this point

Experiment 5: Data Processing 2

5. We also typically calculate four terminal conductance $\frac{dI}{dV} = G_{4T} = \frac{1}{R_{4T}}$ and expressed in units of e^2/h (electron charge squared/Planck constant). Each line in the next graph is data from one file:





6. Calculate the transconductance $\frac{dG_{4T}}{dV_{sg}}$ by numerically differentiating each curve. We typically plot this as an intensity plot vs V_{sg} and B , with color representing the transconductance.



Experiment 6: Next-Level Analysis

7. Use R magic to analyze locations of splittings/crossing/avoided crossings as a function of B and Vsg





Pseudo SCPI Message String for chaining arbitrary number of arguments 2

Example:

setConfig >> Sys >> Time >> hh:mm:ss

setConfig >> Sys >> Date >> yyyy:mm:dd

setSignal >> Magnet >> Target Current >> float (-170 to +1)

Inspiration: SCPI:

<VERB>:<NOUN>:<NOUN>



IEEE Mandated Commands

All SCPI instruments shall implement all the common commands declared mandatory by IEEE 488.2.

Mnemonic	Name	488.2 Section
*CLS	Clear Status Command	10.3
*ESE	Standard Event Status Enable Command	10.10
*ESE?	Standard Event Status Enable Query	10.11
*ESR?	Standard Event Status Register Query	10.12
*IDN?	Identification Query	10.14
*OPC	Operation Complete Command	10.18
*OPC?	Operation Complete Query	10.19
*RST	Reset Command	10.32
*SRE	Service Request Enable Command	10.34
*SRE?	Service Request Enable Query	10.35
*STB?	Read Status Byte Query	10.36
*TST?	Self-Test Query	10.38
*WAI	Wait-to-Continue Command	10.39

SCPI Requirements

IEEE 488.2 describes the syntax of programming and device behavior to a certain level. Further refinement of the syntax and addition rules are imposed by SCPI in order to give instruments a common “look and feel”. The “Syntax and Style” volume of this standard describes the concepts behind SCPI and sets guidelines for originating new commands. A SCPI device shall follow these guidelines and style requirements.

Required Commands

The following commands are required in all SCPI instruments:


Mnemonic	Command Reference Section	Syntax and Style Section
:SYSTem		
:ERRor	21.8	
[:NEXT]?	21.8.3e (see Note 1)	1996
:VERSION?	19.16 (see Note 2)	1991
:STATus	18	5
:OPERation		
[:EVENT]?		
:CONDITION?		
:ENABLE		
:ENABLE?		
:QUESTIONable		
[:EVENT]?		
:CONDITION?		
:ENABLE		
:ENABLE?		
:PRESet		

<https://forums.ni.com/t5/Example-Program-Drafts/Using-Existing-C-Code-or-a-DLL-in-LabVIEW/ta-p/3499233>


<http://www.the-control-freak.com/SCPI/SCPI.htm>

<http://j123b567.github.io/scpi-parser/>

Pseudo SCPI Message String for chaining arbitrary number of arguments 2



Packaging Stuff (2/2-2/3/2019)



Still have error after LV restart and system restart

<https://support.jki.net/hc/en-us/articles/214135643>

Please follow the following VI list and make sure they are closed.

they are may be called by other VIs that are open in LabVIEW.

Name	Expected Path	Actual Path
Logger.Error.lvclass>CreatePrivateEvents.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error>CreatePrivateEvents.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error>CreatePrivateEvents.vi
Logger.Error.lvclass>DestroyPrivateEvents.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error>DestroyPrivateEvents.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error>DestroyPrivateEvents.vi
Logger.Error.lvclass>Error Log Generator Example 1.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error>Error Log Generator Example 1.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error>Error Log Generator Example 1.vi
Logger.Error.lvclass>Log Error.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Log Error.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Log Error.vi
Logger.Error.lvclass>Logger.Error.GetPrivateEvents.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Logger.Error.GetPrivateEvents.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Logger.Error.GetPrivateEvents.vi
Logger.Error.lvclass>Logger.Error.LogError.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Logger.Error.LogError.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Logger.Error.LogError.vi
Logger.Error.lvclass>Logger.Error.TestLauncher.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Logger.Error.TestLauncher.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Logger.Error.TestLauncher.vi
Logger.Error.lvclass>Process.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Process.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\Process.vi
Logger.Error.lvclass>PrivateEvents-Cluster.ctl	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\TypeDefs\PrivateEvents--Cluster.ctl	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\TypeDefs\PrivateEvents--Cluster.ctl
Logger.Error.lvclass>PrivateEvents-Logger.Error.LogError.ctl	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\TypeDefs\PrivateEvents-Logger.Error.LogError.ctl	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger.Error\TypeDefs\PrivateEvents-Logger.Error.LogError.ctl
Logger.lvclass>Logger.TestLauncher.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Logger.TestLauncher.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Logger.TestLauncher.vi
Logger.lvclass>Process.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Process.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Process.vi
Logger.lvclass>Read path.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Read path.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Read path.vi
Logger.lvclass>Write path.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Write path.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\Logger\Write path.vi
LevyLab.lvclass>Handle Error.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\LevyLab\Handle Error.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\LevyLab\Handle Error.vi
LevyLab.lvclass>LevyLab.TestLauncher.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\LevyLab\LevyLab.TestLauncher.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\LevyLab\LevyLab.TestLauncher.vi
LevyLab.lvclass>Process.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\LevyLab\Process.vi	C:\user\patrick\hg\levylab_inst\vsrclSMOs\LevyLab\Process.vi


Troubleshooting Package Building Problems

Sometimes VIPM reports an error during the package building step. Issues that may cause package building to fail and possible solutions are listed below:

- One or more of your source files are missing dependencies. In other words, a VI is calling another VI which is missing. This can happen if your code depends on an add-on toolkit that is not installed with LabVIEW.
 - You have LabVIEW files in your source with duplicate names. All your source filenames should be unique. To work around this issue, you can place VIs in libraries (*.LVLIB files) which provide name-spacing. Having VIs in separate classes also resolves this issue, since classes provide the same name-spacing as *.LVLIBs.
 - You cannot have some of your source files located beneath LabVIEW. In other words, you cannot build your package, if the source is located beneath LabVIEW. VIPM supports [Symbolic Paths](#) so you don't need to locate your source beneath LabVIEW. When VIPM builds your package it replaces your absolute references to the LabVIEW installation location with symbolic paths. This allows your package to be installed under any LabVIEW version installed in any location under Windows.
 - You cannot have some of your source files cross-linking to locations outside your source folder. For example, you have a VI named **test.vi** in your source and you are also calling **test.vi** inside of LabVIEW. This typically happens when you build your package and then install it. Now, when you try to open the source VIs, LabVIEW links to the version which is installed, not your source. - you should try uninstalling your package from LabVIEW before building a new one.

If you need to keep the package installed while you develop your source code, then you should consider adding a prefix to your VIs during build. This can be done automatically by VI Package Builder in the 'Source Files Settings' page. Note, that the installed package VIs must be the renamed version.

- You cannot have VIs open (in memory) while you perform the build.
 - Do not include images, text files or other non labview files inside of LVLIBs or Classes. There have been reported issues with this.
 - If you have several VIs in separate folders and you are trying to build a single package that includes all of those VIs, then you must select the common folder path as your source folder. For example, if you have 'c:\source\component1' and 'c:\source\component2'. Then your source folder must be: 'c:\source'
 - If you are getting Error Code 43 (see image). The problem is that you have a library (lvlib) that has missing VIs inside it. Open the project and navigate to your lvlib project items. Look for instances of VIs that are missing on disk but are still referenced in the lvlib. If you remove these VIs from your library then the package building will complete correctly.




If you are trying to build packages of a class inheritance then there are special considerations:

- When building a package of a child class separate from the parent class (one package for each class) - the child class VIs (in your source) must call the **installed** parent package not the parent class in the source.
 - If you have several classes in separate folders and want to build a single package of all classes, then you must select the common folder path as your source folder. For example, if you have 'c:\source\class1' and 'c:\source\class2'. Then your source folder must be: 'c:\source'

To debug package install issues do the following

- Uninstall all other packages in the the LabVIEW version you are trying to install under. This will provide a clean slate to test your new package.

Switch Build to 2016??



Rename RemoteControl.lvclass:

Rename RemoteControl.lvclass → RemoteControl.v2.lvclass. Get following error:

The file at

'C:\user\patrick\hg\levylab_inst\.lvsr\SMOs\RemoteControl\RemoteControl.v2.lvclass

RemoteControl.v2.ctl' was expected to have the qualified name

'RemoteControl.v2.lvclass:RemoteControl.v2.ctl', but has the qualified name

'RemoteControl.lvclass:RemoteControl.v2.ctl'.

Rename RemoteControl.v2.lvclass → RemoteControl.lvclass. Get following error:

The file at

'C:\user\patrick\hg\levylab_inst\.lvsr\SMOs\RemoteControl\RemoteControl.lvclass\Re

moteControl.ctl' was expected to have the qualified name

'RemoteControl.lvclass:RemoteControl.ctl', but has the qualified name

'RemoteControl.v2.lvclass:RemoteControl.ctl'.



Rename again RemoteControl.lvclass → RemoteControl.lvclass. Seems to work now

Instruments.vipb

****Experimental Build****

[1.1.0.4]

includes:

LevyLab.lvclass

Instrument.lvclass

Configuration.lvclass

Logger.lvclass

Logger.DSC.lvclass

RemoteControl.lvclass

RemoteControl.STM.lvclass

not included:

Logger.syslog.lvclass

RemoteControl.ZMQ.lvclass

Example Instrument.lvclass

Example Instrument UI.lvclass

Instrument's API (override) folder (Go to Angle, etc.)

Dependencies

Joe Albro – 9:08PM on Feb 3

•[Patrick Irvin](#) just in case this helps, all of the dependencies that I needed for the SR7270 Example included:

Remote Repository

•The Instrument Repository

•The G UID Repository

•And I had to install:

Structure Error Handler Package

•DSC (I have this one on other computers but not mine since I reinstalled LabVIEW)

•Syslog Package

•I already had but others may not have:

MGI Error Handler

•JKI JSON

I hope that this may help when it comes to making a package/example project.